# iWard Control Server

DAANTECH	SEMBLE MATES	0 4193 · 0 0053	0 32414 0 1482	O PTAIL	R.
reke 0 -	transfer and the second	113B 2010	1202B 2010	1312D and	1410D (8155)
	17 at 12	EAR.	11 10 10	15.0.15	从共有
0.0.0	(0 2 434 83) C 1207	0 0-2488 C 8352	(I votas) (e sais	0 ******** • ****	(0 8 49 486) (C 8868
	1504A 3180	1605B 3180	1710A (1911	1811A #8#	208B (stat)
	SAT	磺胺碱	HHR	満仕録	680
	(0.010000) (@ 10000	(0 010832)(c 2324	) (0 0 4 4 4 1 (e ####	0 8.40927 (C #308	(if it coutant) (if it to an
	2310B 110	2409C 311	313A 311	3306A 3488	3707A with
	11/64	<b>RAR</b>	610	20,00 12	RRE
	(D. HANGED) (Q. NRAT	(Salasa) (Caste		0	
	38A 311	398 xmt	40C 810	41C 311	42A 9111
	與王卿	15天15月	REA		7,0.9
			(0	(0 8-5937) C 5315	
	43 B 1104	44D 1811	45A 811	4607 minut	478 111

## iWard Control Server datasheet

#### **Features**

- Product classification: 500 clients, 200 clients, 100 clients, 10 clients
- Instant graded demand notification display
- Dynamic management of bed information
- Dynamic adjustment of nursing staff scheduling
- Various statistical reports, historical records query, audit ward performance
- SIP call settings
- System administration management
- SIP Server provide communication function
- Communicate with patient and control center immediately

## Introduction

The iward server is mainly divided into AP server and DB server, which are classified into 500, 200, 100, and 10 according to the number of connections. Provide humanized backend to do bed management, patient management, Nursing schedule, SIP communication settings, upload videos and files. And support Restful API integration The SIP server is classified into 500, 200, 100, and 10 according to the number of connections. Witch provides communication function for client devices.

## **Key Functions**

#### **Nursing schedule**

Nursing scheduling related settings



#### Nursing request list

You can search for nursing notification statistics in the time interval and display the results below



#### System backend

Provide system automatic remote installation, update, delete, upload corresponding logos, and upload files

ADVANTECH		 01111 01111 01		A
		1191		
	DR WARED	 -	and here	
	100			
	a (ac. 12 -			
		AL 010		
	14	31.24	++	++++
		alasia.		
		A REAL		:
		1995		

#### **SIP** management

Provides the function of matching the SIP number of the bed

	RT   ADRES   DRET -				
PRIMER G1000					
00	88.1-98		***	1.0.0	1
ŏ	84-285		+#1		
õ	84-64	94		8517	
00	** ***	1910		-	
00	60				
00	44-44		-	+1	
00	841-88	4.94	204	RSer	
00	841-55		8.8.4	37	
00	84			a2	
00	84**		***	****	
		0			

AD\ANTECH

All product specifications are subject to change without notice.

# **Specifications**


### **Control center**

Service	Function	Description
IDVerification	Login verification	Enter account password or card login to verify user identity, and restrict permission display function.
	Date/Time	Display the current date and time
Information	Nurse station number	Display the nurse station number
Information	logged in account	Display of currently logged in account
	Language	Switch language
	Bed information	Different colors show different genders, male patients are light blue, and female patients are pink. Provide patient name hidden code processing to protect patient privacy. When the patient has a care needs notification item, a reminder pattern is displayed. Show bed number and primary nurse staff.
Bed information	Patient details	Basic information, you can modify patient information such as language and phone For detailed information, you can edit the diagnosis and special notes. The treatment plan shows the patient's treatment date, time, items and status.
	Nursing request	View all requirements notifications
	Inquery function	Filter the list you want to see based on search criteria
	Schedule display	Schedule display
	Schedule inquery	You can select the search conditions by date, and display all or empty beds.
Nursing shift management	Add group	Added information about nursing teams, including day shifts, nursing staff, group, shifts, fire teams, and hospital beds.
	Export function	You can export schedule information for the day or month
	Annousment	Nursing station announcement settings can be set on the control center, and its content can be displayed on the dashboard
	Promotional video	You can upload a video and limit the format of MP4 to provide PIT promotion video playback
	Contact information	Show the contact information of the nursing station medical team and related members of the nursing station
	Device status	Display the device status
Data management	Questionnaire management	Provide management functions for adding, modifying, and deleting satisfaction survey forms
Data manayomont	Education management	Provide management functions for adding, modifying, and deleting health education categories and files
	Living info. management	Provides the classification and file management functions of adding, modifying, and deleting living information
	Environmental Introduction Management	Provide management functions for adding, modifying, and deleting files
	Hospitalization instructions management	Provide management functions for adding, modifying, and deleting files
Nursing request	Nursing request notification	Three notification levels can be displayed (healthcare needs / general needs / other needs).
management	Release notification	Clicking on any of the required notification items will automatically jump to the notified bed, making it easy for caregivers to find information about specific patients.
Nursing Communication	Nursing Communication	Can receive calls from the ward and calls from the nursing app
System management	Account management	Add / modify the account password
Other	Version info.	When you click Logo, the version number and version update time will be displayed in the window.

#### **Backend**

Service	Function	Description
	Parameter	Provide system password (PIT/BIT), preset 0000 Provide SIP communication setting
	Menu	Provide users to add, modify, and delete function menus to determine whether the dashboard or console menu is displayed.
Maintenance management	Communication	Provide users with management functions for adding, modifying, and deleting SIP accounts
(Client)	Communication pairing management	Provide SIP number corresponding bed management function.
	System backend	Provide system automatic remote installation, update, delete, upload corresponding logos, and upload files
	Account management	Add / modify the account password
	Device status	Display the device status
	System managememt	Added and adjusted IWard subsystem.
	Code management	Add, delete, and modify code categories.
Vendor Management (Development side)	Hospital management	Adjust client display settings.
	Module management	Adjust the console and open the background menu function.
	Personnel management	Add, delete, and modify user account and password logins in the background.
	Groug management	Add, delete, and modify nursing group
Language	Language	Switch language

## iWard Architeture



# **Minimum Requirement**

*iWard AP Server		
Hardware	CPU 4 core RAM 8G Up HD C:250G D:250G	
Software	OS: Windows Service 2012R up .Net Framework 4.5.2 Up .Net Core 2.2 Up IIS 8.5+	
*iWard SIP Server		
Hardware	CPU 4 core RAM 4G Up HD C:250G D:250G	
Software	OS: CentOS7+	

*iWard DB Server	
Hardware	CPU 4core RAM 8G Up HD C:250G D:250G
Software	OS: Windows Service 2012R up MSSQL Service 2014 up

# **Ordering Information**

Product Name	Product Number	Description
iWard Control Center	31ASIWSVMSSDE0	iWard Server License for 500 clients
iWard Control Center	31ASIWSVMS20E0	iWard Server License for 200 clients
iWard Control Center	31ASIWSVMS10E0	iWard Server License for 100 clients
iWard Control Center	31ASIWSVMS01E0	iWard Server License for 10 clients

#### \*Option

Product Name	Product Number	Description
Patient Call Service	31ASCMSVLXSDE0	Patient Call SIP Server License for 500 clients
Patient Call Service	31ASCMSVLXSDE0	Patient Call SIP Server License for 200 clients
Patient Call Service	31ASCMSVLXSDE0	Patient Call SIP Server License for 100 clients
Patient Call Service	31ASCMSVLXSDE0	Patient Call SIP Server License for 10 clients